

SOLUTIONS FOR YOUR GOODLIFE!

- 3 layers shoe rack for shoe organizer and storage.
- Rotary shoe drawer.
- High-elastic foam seat cushion with leather cover.
- Knock down package with poly foam and strong carton protecting.
- Easy installing with photo instruction manual.

Design And Color





GLS18815BR



GLS18815BR



GLS18815BR



GLS18815BR



GLS18815BR



GLS18815BR



GLS18815BR

Product Features

Production Name	Wooden color stool shoe cabinet with seat cushion	Brand	Goodlife
Item No.	GLS18815B		
Product size	H42* W63*D30cm (Inch : H16.54*W24.8*D11.81)		
Colors	White, black, brown, wooden color		
Function	floor standing, shoe storage and organizer		
Material	EU-standard MDF with melamine		
Package size	71x46.5x18 cm (Inch:27.95*18.31*7.09)		
Package	KD packing, 1 piece per 1 carton, knock down package, Protected by poly foam, Double layer strong master carton.		
Carton CBM	71*46.5*18 =0.059m ³		N.W./G.W. 14/15kg
20GP	474pcs	40HQ	1152pcs
			MOQ 100pcs

Workmanship And Package

Woodwork

Accessories processing



Constuction

Assembly of sub-pieces
of semi-finished products.



Hand sanded

3 times

Sanding Polishing.



Painting

Automatic spray painting &
Artificial spray painting.



Drop Test

Drop test at International
Standard.



Package

Forming polyfoam ,
postage package.KD package.



Assembly of finished product

MANUAL ASSEMBLY OF SUB-PARTS AND
HARDWARE ACCESSORIES.





Wood Cutting



Drill Holes




Open Slot



sealing with PVC

Material



VS



Particle board

Certificates



Our team





Our Business services:

1. Your inquiry related to our products or prices will be replied in 12hours in working date.
2. Experience sales answer your inquiry and give you related business service.
3. OEM&ODM are welcome, we have over 10 years' experience working with OEM project.
4. We are looking for the sales exclusive of our ODM wooden mirror jewelry cabinet.
5. Sales exclusive agent sales right have been protect.

Package information:

1. Forming polyfoam for the post package
2. 6-side strong polyfoam are used for safe package during transportation.
3. Assembled kit and User Manual are available for each cabinet.

Production lead time

1. 25~50 days after received deposit.
2. Normal on production line production lead time is 15days after received deposit.

Hot Selling



GLD17047



GLD13315



GLD14740



GLD15315



GLD15347



GLD08012



GLD08055



GLD08057



GLD08067





GLS17011WT



GLS17122WH



GLS17017WH



GLS17022WH



GLS18812WT